
## Walker Oxygen Tank Holder JO2W

1. Inventory Oxygen Tank Holder Assembly.







2 Clamps with Hardware



- 2. Remove **Tension Knob** and **Clamps** from their attachment points inside the frame.
- 3. Remove long screw hardware from **Clamps**.
- 4. Attach Clamps to the <u>outside</u> of Cylindrical Frame. There are 5 threaded holes on the back of the Cylindrical Frame. This is to assist in the positioning adjustment of the Cylindrical Frame onto the walker. Using a #2 Phillips screwdriver, <u>Hand tighten</u> the Clamp mounting screws. **DO NOT OVER TIGHTEN**. A thread Lock compound can be used to keep mounting screws from backing out.

NOTE: It is recommended that the Clamps be mounted as far apart as possible. Each application is different so use your best judgment when mounting Clamps onto the Cylindrical Frame.



Completed Assembly

5. Once the Clamps are on the leg of the walker, they can be slid up and down to assist in the positioning adjustment of the Cylindrical Frame where you want it on the walker.

## Walker Oxygen Tank Holder JO2W

- 6. Once the **Cylindrical Frame** is where you want it on the leg of the walker, Re-install long screw hardware into the **Clamps** and tighten down the hardware on the **Clamps**. Hand tight. **DO NOT OVER TIGHTEN**.
- 7. Gently slide an oxygen bottle into the **Cylindrical Frame**.

NOTE: This oxygen tank holder is intended for use with size D or E tanks **ONLY**.

8. Insert the **Tension Knob** into the top of the **Cylindrical Frame** on the <u>outside</u>. <u>Hand tight</u>. **DO NOT OVER TIGHTEN**.



Completed Assembly